

Supplementary Table 1. Oligomer sequences

Oligomer name		Oligomer sequence (5'-3')	Amplicon size
Standard	Telomere standard	(TTAGGG) ¹⁴	84 bp
	36B4 standard	CAGCAAGTGGGAAGGTGTAATCCGCTCCACAGACAAGGCCAGGACTCGTTTG TACCCGTTGATGATAGAATGGG	75 bp
PCR primers	Telomere forward	CGGTTTGTTGGGTTTGGGTTTGGGTTTGGGTTTGGGTT	76 bp
	Telomere reverse	GGCTTGCCTTACCCTTACCCTTACCC TTACCCTTACCCT	
	36B4 forward	CAGCAAGTGGGAAGGTGTAATCC	75 bp
	36B4 reverse	CCCATTCTATCATCAACGGGTACAA	